



FOR IMMEDIATE RELEASE:

CONTACTS:

Matt Smith
Healthmark Industries
(586) 774-7600
msmith@hmark.com

Eliminate Confusion Between Reprocessed and Soiled Devices

Fraser, MI, January 17, 2016 - Healthmark has announced SST Liners to their product line. These colored coded liners help to eliminate confusion between reprocessed and soiled medical devices.

The single-use liners assist with identification, temporary storage and transportation of medical equipment by providing an enclosed environment for instruments. Simply, line the base of the tray using the clear liner to ensure interior base and walls are covered; place clean device in tray along with the red biohazard liner for post-procedure use and cover with the green clean liner. Upon procedure completion discard the green liner, place used device back in tray and cover with the red biohazard liner to indicate item is ready for transport to the Reprocessing area.

Made from low-density polyethylene, the color-coded tray liners are 36 x 18 Inches, 1 Mil and available in the following colors: The SST Liners are available in three colors and available in a box of one-hundred.

About Healthmark Industries Co., Inc.

Since 1969, Healthmark Industries Company, Inc. has developed and marketed innovative solutions to aid healthcare facilities in their delivery of surgical instruments and other life saving medical devices to patients. Healthmark Industries mission is to continue to innovate, continue to support and continue to serve the healthcare provider industry and support services that make it possible to deliver quality healthcare. Visit www.hmark.com for more information.